



PREempt™ CS20 USE & TESTING RESULTS LOG

PREempt CS20 is good for reuse in a soaking tray for up to 14 days.

Test CS20 solution for efficacy on day 7 of use and whenever there is visible debris in the soaker.

ONLY use AHP7 Indicator Test Strips to test CS20. Record the date of testing and the test results in the "Debris Test Date & Result" column.

Please include your initials when entering dates.

Fresh Solution Date	Replace Solution Date (14 days)	Day 7 Test Date & Result	Debris Test Date & Result	Debris Test Date & Result

Fresh Solution Date	Replace Solution Date (14 days)	Day 7 Test Date & Result	Debris Test Date & Result	Debris Test Date & Result