


TheraMax

JE-7000



TheraMax

Thermo acupressure massager

TheraMax -Thermo Acupressure Massager

TheraMax is solely committed to help and fulfill each TheraMax user's desire to be healthy. TheraMax is a unique acupressure system engineered by fine and recognized doctors, engineers, and mechanics. This therapeutic massage bed with the superbly engineered Twin Jades is a revolution to Therapeutic Massage Bed industry.

The newly created and the first ever introduced sharp premium twin jades provide much better pressure than any other competitors' devices. Heat is essential key to having a successful therapy and with the newly created twin premium jades, TheraMax users can receive more heat which creates much better result.



Therapeutic Effect & Efficacy

With the new concept of thermo acupressure massager, TheraMax is designed by combining Oriental Medical principle such as moxibustion and acupressure with Western medical technology, the Chiropractic.

Chiropractic and acupressure effect

TheraMax Internal Jades that are located in upper mat travel the spine line from cervical region to sacrum (hips). It helps to maintain and improve alignments and increase flexibility of entire spine. The jades stimulate the vital acupuncture points to boost metabolism while traveling through the spine. It also circulates the body energy, which improves immune system.

Heat Therapy (Moxibustion function) Effect

When the heat penetrates and relaxes the muscle tissues, it may alleviate the pain and fatigue. Heat therapy stimulates the capillary vessel to extend and improve blood stream.

Massage Effect

TheraMax premium rolling jade is superior to any other devices in relieving myalgia (muscle pain). The internal jade automatically massages muscles and tendons around the spine while traveling through the spine. By massaging the upper and lower parts of body area, it relaxes hardened nerve roots and improves the circulation of energy.


Far Infrared Ray from Jade & Bio-Ceramic

Infrared Ray is one of the heat rays, which are invisible and have a longer wavelength than visible red heat rays. TheraMax infrared radiant energy directly penetrates the body to a depth of over 1.5 inches. Far Infrared Ray also provides incredible benefits to your body. It stimulates the cell function, discharge accumulated toxins, increase blood circulation, and loosen all tight muscles. For more than 30 years, far infrared rays have been beneficially used to treat patients. The use of these rays for medical purpose has proven to be effective in many physical ailments.

Comfortable and soft ergonomically curved design

Ergonomics Design
TheraMax's lower mat applies Leg Guard System. It keeps your legs straight on the lower mat. It is applied to the mat to support human body in comfortable and pleasant condition with a beautiful appearance and high quality materials.


Diagnostic Control System
TheraMax was designed with a diagnostic control circuit, so you can find technical problem or any malfunction from a display-monitor that shows accurate messages.



SPECIFICATION

| | |
|-----------------------|----------------|
| Voltage | AC 120V / 230V |
| Watt | 130 W Maxium |
| Frequency | 60 Hz |
| Operating Temperature | 45-70 C degree |

- 1 Internal Heater**
It moves from head to hips giving thermal and acupressure effect. The internal Projector, with it's heated jade massage rollers, moves slowly along the spinal contour providing massage, acupressure, and far-infrared light and heat therapy to the back. Major acupuncture meridians are stimulated, energizing and revitalizing the body.
- 2 6 Ball Projector**
It is used to deliver the effects of far infrared rays and acupressure on the stomach. The ergonomically designed jade, 6-ball projector provides far-infrared heat and light therapy and may be conveniently used on any part of the body.
- 3 Bio Ceramic**
This special, auxiliary mat of Bio Ceramic composition emits healthy negative ions, and is used for enhanced far-infrared heat and light therapy.
- 4 Safe Connector**
Connector was designed to prevent any potential accident.
- 5 Remote Controller**
All the functional settings are easily adjusted by the remote control.
- 6 Wooden Frame**
Classical frame is designed elegance shapes and made with wood material solidly.
- 2 Ball Projector (optional)**
It is used to deliver the effects of far infrared rays and acupressure on the neck, face.
- Leg Roller (optional)**
The leg massage roller is designed by human mechanism, ergonomics. It provides far-infrared ray and heat therapy with soft internal jades' movement.
- High Technology**
TheraMax JE-7000 was developed with famous scientists and engineers, have developed the world's best thermal massager.



Authorized Supplier
Fernanda's Beauty Products
1225 Aerowood Drive
Mississauga, Ontario,
Canada L4W 1B9
Toll Free: 1-800-862-1447
www.fernandas.com